



CYSTOSCOPY



FULL HD CYSTOSCOPES

- Colour coded rings identify the direction of view
- 3-step light post adapters for the most common light cables
- Sapphire glasses to protect lens system
- Optimized glass fibre arrangements for an uniform image brightness
- Fully autoclavable at 134°C and 2.3bar



FULL HD

Ø 4mm Endoscopes for Cystoscopy



● 0° DIRECTION OF VIEW

300mm working length
356mm overall length 253-804-000*



○ 12° DIRECTION OF VIEW

300mm working length
356mm overall length 253-804-012*



● 30° DIRECTION OF VIEW


300mm working length
356mm overall length 253-804-030*



● 70° DIRECTION OF VIEW


300mm working length
356mm overall length 253-804-070*

Sheaths for use with Ø4mm Cystoscopes




Standard Cysto-Urethroscope sheaths with 2 lateral stopcocks and standard obturator with atraumatic tip

- 253-300-017* Cysto-Urethroscope Sheath 17 Charr. ●
- 253-300-019* Cysto-Urethroscope Sheath 19 Charr. ●
- 253-300-020* Cysto-Urethroscope Sheath 20 Charr. ●
- 253-300-022* Cysto-Urethroscope Sheath 22 Charr. ●
- 253-300-025* Cysto-Urethroscope Sheath 25 Charr. ○



Long Window Cysto-Urethroscope sheaths with 2 lateral stopcocks and standard obturator with atraumatic tip

- 253-300-119 Cysto-Urethroscope Sheath 19 Charr. ●
- 253-300-120 Cysto-Urethroscope Sheath 20 Charr. ●
- 253-300-122 Cysto-Urethroscope Sheath 22 Charr. ●
- 253-300-125 Cysto-Urethroscope Sheath 25 Charr. ○



Central Valve Cysto-Urethroscope sheaths with 2 stopcocks and central valve
Sheaths with standard obturator and atraumatic tip

- 253-310-017 Cysto-Urethroscope Sheath 17 Charr. ●
- 253-310-019 Cysto-Urethroscope Sheath 19 Charr. ●
- 253-310-020 Cysto-Urethroscope Sheath 20 Charr. ●
- 253-310-022 Cysto-Urethroscope Sheath 22 Charr. ●
- 253-310-025 Cysto-Urethroscope Sheath 25 Charr. ○

Central valve for easy regulation
of suction and irrigation



Obturator for use with Cystoscope sheaths



Standard Obturators with atraumatic tip and locking mechanism for all RZ Cysto-Urethroscope sheaths

- 253-000-117 Obturator 17 Charr. ●
- 253-000-119 Obturator 19 Charr. ●
- 253-000-120 Obturator 21 Charr. ●
- 253-000-122 Obturator 23 Charr. ●
- 253-000-125 Obturator 25 Charr. ○




Visual Obturators with locking mechanism for all 4mm RZ Cystoscopes

- 253-000-217 Obturator 17 Charr. ●
- 253-000-219 Obturator 19 Charr. ●
- 253-000-220 Obturator 20 Charr. ●
- 253-000-222 Obturator 22 Charr. ●
- 253-000-225 Obturator 25 Charr. ○




253-315-000
Hamada Cystoscope Adapter for use with RZ Cystoscope sheath

Cystoscope Bridge for use with Cysto-Urethroscope sheaths




- 1** 253-000-100* without channel
- 2** 253-000-101* with 1 channel for 1 instrument with max. 12 Charr.
- 3** 253-000-102* with 2 channels for 1 instrument with max. 12 Charr. or 2 instruments with max. 8 Charr.

Albarran deflecting mechanism for use with 19-25 Charr. RZ Cysto-Urethroscope sheaths



- 253-000-105* Albarran deflecting mechanism with 2 working channels **and ratchet**
- 253-000-107 Albarran deflecting mechanism with 1 working channel **and ratchet** (without illustration)




- 253-000-106 Albarran deflecting mechanism with 2 working channels **without ratchet** (without illustration)
- 253-000-108 Albarran deflecting mechanism with 1 working channel **without ratchet**








Sheath - Instrument combinations

SHEATHS	COLOR	ALBARRAN		BRIDGE	
		1 working channel	2 working channels	1 working channel	2 working channels
25 Charr.	○	1 x 10 Charr.	2 x 8 Charr.	1 x 12 Charr.	2 x 8 Charr.
22 Charr.	●	1 x 9 Charr.	2 x 6 Charr.	1 x 10 Charr.	2 x 7 Charr.
20 Charr.	●	1 x 6 Charr.	2 x 5 Charr.	1 x 7 Charr.	2 x 6 Charr.
19 Charr.	●	1 x 5 Charr.	2 x 4 Charr.	1 x 6 Charr.	2 x 5 Charr.
17 Charr.	●	-	-	1 x 5 Charr.	-

Flexible Instruments for use with cystoscopes





Working length: 400mm
Rotatable and detachable


	Flexible	Semi-rigid
 Grasping Forceps, serrated jaws, double action	351-114-105* 5 Charr. 351-114-107* 7 Charr. 351-114-109 9 Charr.	351-114-305 5 Charr. 351-114-307 7 Charr.
 Biopsy Forceps, oval cup, double action	351-120-105* 5 Charr. 351-120-107* 7 Charr. 351-120-109 9 Charr.	351-120-305 5 Charr. 351-120-307 7 Charr.
 Serrated Biopsy Forceps, oval cup with aggressive serration, double action	351-124-107 7 Charr.	351-124-307 7 Charr.
 Scissors, single action	351-100-105* 5 Charr. 351-100-107 7 Charr. 351-100-109 9 Charr.	351-100-305 5 Charr. 351-100-307 7 Charr.
 DJ Forceps, double action	351-107-107 7 Charr. 351-107-109 9 Charr.	351-107-307 7 Charr.
 Grasping Forceps, serrated jaw, double action, 730mm working length	253-958-105 5 Charr.	
 Biopsy Forceps, double action, 730mm working length	253-957-105 5 Charr.	

OPTICAL FORCEPS

- for use with \varnothing 4mm Cystoscopes



Scope is **below** forceps



Scope is **above** Forceps



253-420-100*
Spoon Biopsy Forceps with round cup,
double action



253-421-100
Grasping Forceps with serrated jaws,
double action



253-422-100
Biopsy Punch Forceps, thru-cutting,
double action



253-423-100
Scissors, hook type blades,
double action



253-426-100*
Mini Biopsy forceps with oval cup,
double action
for use with 0° cystoscope




253-431-100*
Biopsy forceps with **large** oval cup,
double action
for use with 30° cystoscope



253-425-100
Biopsy forceps with mini oval cup,
double action
for use with 4mm, 70° cystoscope




253-427-100*
Stent Grasping forceps
with serrated jaws, double action
for removal of stents




253-429-100* Adapter for use of optical forceps with resectoscope sheaths


FLEXIBLE INSTRUMENTS



253-940-005 3-Prongs Grasper, flexible, 600mm working length, 5 Charr.




253-585-005 Dormia Metal Stone Basket, flexible, 4 wires, 600mm working length, 5 Charr.




Dormia Nitinol Stone Basket, 600mm working length, 5 Charr.

- 253-581-045 with 4 wires, complete with handle
- 253-582-045 Spare Basket only, 4 wires
- 253-582-055 Spare Basket only, 5 wires
- 253-582-065 Spare Basket only, 6 wires




Polyp Snare with teflon tube, 500mm working length with HF connection, 7 Charr.

253-580-007 with handle, complete
 253-581-007 Spare Snare only




253-584-005 Zeiss Loop Catheter, package with 5 pieces, sterile, single use, 5 Charr.




253-000-235 Injection Cannula, semi-rigid, 4 Charr.
 distal end 2 Charr., with rotatable Luer-Lock adapter




FLEXIBLE ELECTRODES AND HF CABLES

Ball Electrode, 530mm working length




	253-370-003	3 Charr.
	253-370-103	3 Charr., 115cm working length
	253-370-005*	5 Charr.
	253-370-007	7 Charr.
	253-370-010	10 Charr.



Needle Electrode, 530mm working length

	253-371-003	3 Charr.
	253-371-005	5 Charr.
	253-371-007	7 Charr.


Loop Electrode 530mm working length

	253-372-005	5 Charr.
	253-372-007	7 Charr.
	253-372-010	10 Charr.



Monopolar High Frequency cables, 4,5mm length


253-610-001*	with 5 mm plug for Erbe HF-unit
253-610-002*	with 4 mm plug for Martin HF-unit
253-610-003*	with 8 mm plug for Valleylab HF-unit
253-610-004*	with 4 mm plug for older HF-units
700-501-453*	SDS plug



CYSTOSCOPY ACCESSORIES



301-701-100* Grey rubber cap for use on all working channels



253-990-800 French Catheter Scale



Cleaning Brush, flexible

253-961-103 \varnothing 3mm, 100mm length,
for use in 2.0 - 2.6mm working channels

253-961-102 \varnothing 2.5mm, 120mm length,
for 5 Charr. working channels



300-215-000

Cleaning Brush, \varnothing 13mm, 60mm brush length, 230mm overall length