


Cleaning and Sterilization



STERILIZATION CONTAINER

CONTAINER MODELS



STANDARD

model

- Easy insertion and change of filter materials and filter retainer
- Extra long handles with silicone
- Stable and safe closing mechanism canceling overpressure
- Tag holders for clear traceability



SAFETY

model

- Easy insertion and change of filter materials and filter retainer
- Extra long handles with silicone
- Stable and safe closing mechanism canceling overpressure
- Tag holders for clear traceability
- Additional safety lid with exterior press-button



ADVANCED SAFETY

model

- Easy insertion and change of filter materials and filter retainer
- Extra long handles with silicone
- Stable and safe closing mechanism canceling overpressure
- Tag holders for clear traceability
- Color matched lid and handles
- Integrated safety lid
- Increased stability due to reinforced container walls
- Integrated spring system

SYMBOL EXPLANATION



Bottom: non-perforated



Bottom: perforated



Safety lid
Bottom: non-perforated



Safety lid
Bottom: perforated



Bottom only: non-perforated



Bottom only: perforated



Standard lid



Safety lid



Advanced Safety lid

LID COLORS



● Silver
End of numeral: 1



● Blue
End of numeral: 2



● Red
End of numeral: 3



● Green
End of numeral: 4



● Yellow
End of numeral: 5



● Black
End of numeral: 6

STERILIZATION CONTAINER

Suitability of our sterile containers for steam sterilization procedures¹⁾ as per DIN 58953 part 9 (issue Mai 1987)

	Packing		Steam sterilization procedure acc. to DIN 58946, part 1				
	Description	Symbolic representation and packing system	Gravitation procedure	Pre-vacuum procedure	Fractional vacuum procedure	Fractional flow procedure	Steam injection procedure
with filters (StB-F)	Textiles in sterilization container with filter in lid and bottom acc. to DIN 58952 part 1		✓	✓	✓	✓	✓
			✓	✓	✓	✓	✓
	Textiles in sterilization container with filter in lid, bottom not perforated acc. to DIN 58952 part 1		✗	✗	✓	✓	✓
			✗	✗	✓	✓	✓
	Textiles in sterilization container with filter bottom, lid not perforated acc. to DIN 58952 part 1		✓	✓	✓	✓	✓
			✗	✓	✓	✓	✓
	Instruments in wire basket acc. to DIN 58952 part III, in sterilization container with filter in lid and bottom acc. to DIN 58952 part 1		✓	✓	✓	✓	✓
	Instruments in wire basket acc. to DIN 58952 part III, in sterilization container with filter in lid, bottom not perforated acc. to DIN 58952 part 1		✗	✓	✓	✓	✓
	Instruments in wire basket acc. to DIN 58952 part III, in sterilization container with filter in bottom acc. to DIN 58952 part 1		✓	✓	✓	✓	✓
	with valves (StB-V)	Textiles in sterilization container with valves in lid and bottom acc. to DIN 58952 part 1		✗	✓	✓	✓
			✗	✗	✓	✓	✓
Textiles in sterilization container with valves in lid, no valves on the bottom			✗	✗	✓	✓	✗
			✗	✗	✓	✓	✗
Instruments in wire basket in sterilization container with valve in lid and bottom			✗	✓	✓	✓	✓
Instruments in wire basket in sterilization container with valve in lid, bottom not perforated			✗	✓	✓	✓	✓

¹⁾ In order to avoid mistakes with regards to service, operation and/or when loading the sterilizer, please refer to table 1 for correct or wrong combination of sterilization containers for standardized steam sterilization procedure.

✓ Applicable ✓ Limited applicable ✗ not applicable

²⁾ Suitable are only sterilization containers guaranteeing evaporation or discharge of condensate.

³⁾ Gravitation procedures having no drying device do not guarantee a sufficient after-drying effect.


1/1 CONTAINER

580 x 280mm



STERILIZATION CONTAINER

1/1 ADVANCED SAFETY CONTAINER



Size (LxWxH)



580x280x100 mm

▲ 100mm

●	600-580-101	600-579-101
●	600-580-102	600-579-102
●	600-580-103	600-579-103
●	600-580-104	600-579-104
●	600-580-105	600-579-105
●	600-580-106	600-579-106



580x280x135 mm

▲ 135mm

●	600-580-131	600-579-131
●	600-580-132	600-579-132
●	600-580-133	600-579-133
●	600-580-134	600-579-134
●	600-580-135	600-579-135
●	600-580-136	600-579-136



580x280x150 mm

▲ 150mm

●	600-580-151	600-579-151
●	600-580-152	600-579-152
●	600-580-153	600-579-153
●	600-580-154	600-579-154
●	600-580-155	600-579-155
●	600-580-156	600-579-156



580x280x200 mm

▲ 200mm

●	600-580-201	600-579-201
●	600-580-202	600-579-202
●	600-580-203	600-579-203
●	600-580-204	600-579-204
●	600-580-205	600-579-205
●	600-580-206	600-579-206



580x280x260 mm

▲ 260mm

●	600-580-261	600-579-261
●	600-580-262	600-579-262
●	600-580-263	600-579-263
●	600-580-264	600-579-264
●	600-580-265	600-579-265
●	600-580-266	600-579-266

1/1 SAFETY CONTAINER



Size (LxWxH)



580x280x100 mm	● 600-583-101	600-584-101
	● 600-583-102	600-584-102
	● 600-583-103	600-584-103
	● 600-583-104	600-584-104
	● 600-583-105	600-584-105
▲ 100mm	● 600-583-106	600-584-106



580x280x135 mm	● 600-583-131	600-584-131
	● 600-583-132	600-584-132
	● 600-583-133	600-584-133
	● 600-583-134	600-584-134
	● 600-583-135	600-584-135
▲ 135mm	● 600-583-136	600-584-136



580x280x150 mm	● 600-583-151	600-584-151
	● 600-583-152	600-584-152
	● 600-583-153	600-584-153
	● 600-583-154	600-584-154
	● 600-583-155	600-584-155
▲ 150mm	● 600-583-156	600-584-156



580x280x200 mm	● 600-583-201	600-584-201
	● 600-583-202	600-584-202
	● 600-583-203	600-584-203
	● 600-583-204	600-584-204
	● 600-583-205	600-584-205
▲ 200mm	● 600-583-206	600-584-206



580x280x260 mm	● 600-583-261	600-584-261
	● 600-583-262	600-584-262
	● 600-583-263	600-584-263
	● 600-583-264	600-584-264
	● 600-583-265	600-584-265
▲ 260mm	● 600-583-266	600-584-266

STERILIZATION CONTAINER

1/1 STANDARD CONTAINER

Size (LxWxH)



580x280x100 mm

▲ **100mm**

●	600-581-101	600-582-101
●	600-581-102	600-582-102
●	600-581-103	600-582-103
●	600-581-104	600-582-104
●	600-581-105	600-582-105
●	600-581-106	600-582-106



580x280x135 mm

▲ **135mm**

●	600-581-131	600-582-131
●	600-581-132	600-582-132
●	600-581-133	600-582-133
●	600-581-134	600-582-134
●	600-581-135	600-582-135
●	600-581-136	600-582-136



580x280x150 mm

▲ **150mm**

●	600-581-151	600-582-151
●	600-581-152	600-582-152
●	600-581-153	600-582-153
●	600-581-154	600-582-154
●	600-581-155	600-582-155
●	600-581-156	600-582-156



580x280x200 mm

▲ **200mm**

●	600-581-201	600-582-201
●	600-581-202	600-582-202
●	600-581-203	600-582-203
●	600-581-204	600-582-204
●	600-581-205	600-582-205
●	600-581-206	600-582-206



580x280x260 mm

▲ **260mm**

●	600-581-261	600-582-261
●	600-581-262	600-582-262
●	600-581-263	600-582-263
●	600-581-264	600-582-264
●	600-581-265	600-582-265
●	600-581-266	600-582-266

1/1



Bottom perforated

Size (LxWxH)



580x280x100 mm	600-585-100
580x280x135 mm	600-585-135
580x280x150 mm	600-585-150
580x280x200 mm	600-585-200
580x280x260 mm	600-585-260



Bottom non-perforated

Size (LxWxH)




580x280x100 mm	600-586-100
580x280x135 mm	600-586-135
580x280x150 mm	600-586-150
580x280x200 mm	600-586-200
580x280x260 mm	600-586-260



Standard Lid

Size (LxW)



580x280 mm	600-587-281
	600-587-282
	600-587-283
	600-587-284
	600-587-285
	600-587-286



Safety Lid


Size (LxW)



580x280 mm	600-588-281
	600-588-282
	600-588-283
	600-588-284
	600-588-285
	600-588-286

STERILIZATION TRAYS

1/1




Screen Basket
with Lid

Size (LxWxH)

545x255x35mm	600-540-535
545x255x55mm	600-540-555
545x255x75mm	600-540-575
545x255x105mm	600-540-605

485x255x35mm	600-485-535
485x255x55mm	600-485-555
485x255x75mm	600-485-575
485x255x105mm	600-485-605



Lid only

545x255mm	600-540-500
485x255mm	600-485-500



Bottom only

545x255x30mm	600-540-530
545x255x50mm	600-540-550
545x255x70mm	600-540-570
545x255x100mm	600-540-600

485x255x30mm	600-485-530
485x255x50mm	600-485-550
485x255x70mm	600-485-570
485x255x100mm	600-485-600


1/1



Rippled wire mesh


Size (LxWxH)

485x255x30 mm	600-485-030
485x255x50 mm	600-485-050
485x255x70 mm	600-485-070
485x255x100 mm	600-485-100
540x255x30 mm	600-540-030
540x255x50 mm	600-540-050
540x255x70 mm	600-540-070
540x255x100 mm	600-540-100



Wire Divider

One Partition	510x250x130 mm	600-901-130
	510x250x180 mm	600-901-180
	510x250x230 mm	600-901-230
Two Partitions	510x250x130 mm	600-901-131
	510x250x180 mm	600-901-181
	510x250x230 mm	600-901-231



Drape Retainer

545x255mm

600-900-003


3/4 CONTAINER

465 x 280mm



3/4 ADVANCED SAFETY CONTAINER

Size (LxWxH)



465x280x100 mm

●	600-469-101	600-470-101
●	600-469-102	600-470-102
●	600-469-103	600-470-103
●	600-469-104	600-470-104
●	600-469-105	600-470-105
●	600-469-106	600-470-106

▲ **100mm**



465x280x135 mm

●	600-469-131	600-470-131
●	600-469-132	600-470-132
●	600-469-133	600-470-133
●	600-469-134	600-470-134
●	600-469-135	600-470-135
●	600-469-136	600-470-136

▲ **135mm**



465x280x150 mm

●	600-469-151	600-470-151
●	600-469-152	600-470-152
●	600-469-153	600-470-153
●	600-469-154	600-470-154
●	600-469-155	600-470-155
●	600-469-156	600-470-156

▲ **150mm**

STERILIZATION CONTAINER

3/4 SAFETY CONTAINER



Size (LxWxH)



465x280x100 mm

▲ 100mm

●	600-463-101	600-464-101
●	600-463-102	600-464-102
●	600-463-103	600-464-103
●	600-463-104	600-464-104
●	600-463-105	600-464-105
●	600-463-106	600-464-106



465x280x135 mm

▲ 135mm

●	600-463-131	600-464-131
●	600-463-132	600-464-132
●	600-463-133	600-464-133
●	600-463-134	600-464-134
●	600-463-135	600-464-135
●	600-463-136	600-464-136



465x280x150 mm

▲ 150mm

●	600-463-151	600-464-151
●	600-463-152	600-464-152
●	600-463-153	600-464-153
●	600-463-154	600-464-154
●	600-463-155	600-464-155
●	600-463-156	600-464-156

3/4 STANDARD CONTAINER



Size (LxWxH)



465x280x100 mm

▲ **100mm**

●	600-461-101	600-462-101
●	600-461-102	600-462-102
●	600-461-103	600-462-103
●	600-461-104	600-462-104
●	600-461-105	600-462-105
●	600-461-106	600-462-106



465x280x135 mm

▲ **135mm**

●	600-461-131	600-462-131
●	600-461-132	600-462-132
●	600-461-133	600-462-133
●	600-461-134	600-462-134
●	600-461-135	600-462-135
●	600-461-136	600-462-136



465x280x150 mm

▲ **150mm**

●	600-461-151	600-462-151
●	600-461-152	600-462-152
●	600-461-153	600-462-153
●	600-461-154	600-462-154
●	600-461-155	600-462-155
●	600-461-156	600-462-156

STERILIZATION CONTAINER & TRAYS

3/4



Size (LxWxH)




465x280x100 mm	600-465-100
465x280x135 mm	600-465-135
465x280x150 mm	600-465-150




Size (LxWxH)



465x280x100 mm	600-466-100
465x280x135 mm	600-466-135
465x280x150 mm	600-466-150



Size (LxW)



465x280 mm	<ul style="list-style-type: none"> ● 600-467-281 ● 600-467-282 ● 600-467-283 ● 600-467-284 ● 600-467-285 ● 600-467-286
------------	---




Size (LxW)



465x280 mm	<ul style="list-style-type: none"> ● 600-468-281 ● 600-468-282 ● 600-468-283 ● 600-468-284 ● 600-468-285 ● 600-468-286
------------	---


3/4



Screen Basket
with Lid


Size (LxWxH)

405x255x35mm	600-405-535
405x255x55mm	600-405-555
405x255x75mm	600-405-575
405x255x105mm	600-405-605




Lid only

405x255mm	600-405-500
-----------	-------------




Bottom only

405x255x30mm	600-405-530
405x255x50mm	600-405-550
405x255x70mm	600-405-570
405x255x100mm	600-405-600



Rippled
Wire Mesh

405x255x30mm	600-405-030
405x255x50mm	600-405-050
405x255x70mm	600-405-070
405x255x100mm	600-405-100



Drape Retainer

405x255mm	600-463-200
-----------	-------------


1/2 CONTAINER

285 x 280mm



1/2 ADVANCED SAFETY CONTAINER

Size (LxWxH)



285x280x100 mm



600-280-101

600-279-101



600-280-102

600-279-102



600-280-103

600-279-103



600-280-104

600-279-104



600-280-105

600-279-105

▲ 100mm



600-280-106

600-279-106



285x280x135 mm



600-280-131

600-279-131



600-280-132

600-279-132



600-280-133

600-279-133



600-280-134

600-279-134



600-280-135

600-279-135

▲ 135mm



600-280-136

600-279-136



285x280x150 mm



600-280-151

600-279-151



600-280-152

600-279-152



600-280-153

600-279-153



600-280-154

600-279-154



600-280-155

600-279-155

▲ 150mm



600-280-156

600-279-156



285x280x200 mm



600-280-201

600-279-201



600-280-202

600-279-202



600-280-203

600-279-203



600-280-204

600-279-204



600-280-205

600-279-205

▲ 200mm



600-280-206

600-279-206



285x280x260 mm



600-280-261

600-279-261



600-280-262

600-279-262



600-280-263

600-279-263



600-280-264

600-279-264



600-280-265

600-279-265

▲ 260mm



600-280-266

600-279-266

STERILIZATION CONTAINER

1/2 SAFETY CONTAINER



Size (LxWxH)



285x280x100 mm

▲ 100mm

●	600-283-101	600-284-101
●	600-283-102	600-284-102
●	600-283-103	600-284-103
●	600-283-104	600-284-104
●	600-283-105	600-284-105
●	600-283-106	600-284-106



285x280x135 mm

▲ 135mm

●	600-283-131	600-284-131
●	600-283-132	600-284-132
●	600-283-133	600-284-133
●	600-283-134	600-284-134
●	600-283-135	600-284-135
●	600-283-136	600-284-136



285x280x150 mm

▲ 150mm

●	600-283-151	600-284-151
●	600-283-152	600-284-152
●	600-283-153	600-284-153
●	600-283-154	600-284-154
●	600-283-155	600-284-155
●	600-283-156	600-284-156



285x280x200 mm

▲ 200mm

●	600-283-201	600-284-201
●	600-283-202	600-284-202
●	600-283-203	600-284-203
●	600-283-204	600-284-204
●	600-283-205	600-284-205
●	600-283-206	600-284-206





285x280x260 mm

▲ 260mm

●	600-283-261	600-284-261
●	600-283-262	600-284-262
●	600-283-263	600-284-263
●	600-283-264	600-284-264
●	600-283-265	600-284-265
●	600-283-266	600-284-266

1/2 STANDARD CONTAINER

Size (LxWxH)



285x280x100 mm

▲ 100mm

●	600-281-101	600-282-101
●	600-281-102	600-282-102
●	600-281-103	600-282-103
●	600-281-104	600-282-104
●	600-281-105	600-282-105
●	600-281-106	600-282-106



285x280x135 mm

▲ 135mm

●	600-281-131	600-282-131
●	600-281-132	600-282-132
●	600-281-133	600-282-133
●	600-281-134	600-282-134
●	600-281-135	600-282-135
●	600-281-136	600-282-136



285x280x150 mm

▲ 150mm

●	600-281-151	600-282-151
●	600-281-152	600-282-152
●	600-281-153	600-282-153
●	600-281-154	600-282-154
●	600-281-155	600-282-155
●	600-281-156	600-282-156



285x280x200 mm

▲ 200mm

●	600-281-201	600-282-201
●	600-281-202	600-282-202
●	600-281-203	600-282-203
●	600-281-204	600-282-204
●	600-281-205	600-282-205
●	600-281-206	600-282-206



285x280x260 mm


▲ 260mm

●	600-281-261	600-282-261
●	600-281-262	600-282-262
●	600-281-263	600-282-263
●	600-281-264	600-282-264
●	600-281-265	600-282-265
●	600-281-266	600-282-266

STERILIZATION CONTAINER


1/2

Screen Basket
with Lid




Size (LxWxH)

255x245x35mm	600-255-535
255x245x55mm	600-255-555
255x245x75mm	600-255-575
255x245x105mm	600-255-605



Lid only

255x245mm	600-255-500
-----------	-------------



Bottom only

255x245x30mm	600-255-530
255x245x50mm	600-255-550
255x245x70mm	600-255-570
255x245x100mm	600-255-600

1/2


Rippled
Wire Mesh



Size (LxWxH)

255x245x30mm	600-255-030
255x245x50mm	600-255-050
255x245x70mm	600-255-070
255x245x100mm	600-255-100

Wire
Divider



255x245x130mm	600-901-132
255x245x180mm	600-901-182
255x245x280mm	600-901-232

Drape Retainer



255x255mm

600-288-300

SPECIAL SIZE CONTAINER



Size (LxWxH)



600x400x120 mm

●	600-600-121	600-601-121
●	600-600-122	600-601-122
●	600-600-123	600-601-123
●	600-600-124	600-601-124
●	600-600-125	600-601-125
●	600-600-126	600-601-126

▲ 120mm



600x400x180 mm

●	600-600-181	600-601-181
●	600-600-182	600-601-182
●	600-600-183	600-601-183
●	600-600-184	600-601-184
●	600-600-185	600-601-185
●	600-600-186	600-601-186

▲ 180mm



470x337x51 mm
533x337x101 mm

600-470-050
600-533-100

Square Hole
Wire Mesh

STERILIZATION CONTAINER

Size (LxWxH)



285x280x55 mm

▲ 55mm

●	600-271-051
●	600-271-052
●	600-271-053
●	600-271-054
●	600-271-055
●	600-271-056



▲ 85mm

●	600-271-081	600-272-081
●	600-271-082	600-272-082
●	600-271-083	600-272-083
●	600-271-084	600-272-084
●	600-271-085	600-272-085
●	600-271-086	600-272-086


Square hole
wire mesh




255x245x30mm	600-255-030
255x245x50mm	600-255-050

Size (LxWxH)

Screen Basket
with Lid



255x245x35mm	600-255-535
255x245x55mm	600-255-555
255x245x75mm	600-255-575
255x245x105mm	600-255-605




Lid only

255x245mm	600-255-500
-----------	-------------



Bottom only

255x245x30mm	600-255-530
255x245x50mm	600-255-550
255x245x70mm	600-255-570
255x245x100mm	600-255-600



Drape Retainer

255x255mm	600-288-300
-----------	-------------

STERILIZATION CONTAINER



Size (LxWxH)



300x140x40 mm

●	600-141-041	600-142-041
●	600-141-042	600-142-042
●	600-141-043	600-142-043
●	600-141-044	600-142-044
●	600-141-045	600-142-045
●	600-141-046	600-142-046

▲ 40mm



300x140x70 mm

●	600-141-071	600-142-071
●	600-141-072	600-142-072
●	600-141-073	600-142-073
●	600-141-074	600-142-074
●	600-141-075	600-142-075
●	600-141-076	600-142-076

▲ 70mm



300x140x100 mm

●	600-141-101	600-142-101
●	600-141-102	600-142-102
●	600-141-103	600-142-103
●	600-141-104	600-142-104
●	600-141-105	600-142-105
●	600-141-106	600-142-106

▲ 100mm



Size (LxWxH)



300x140x40 mm	600-143-040
300x140x70 mm	600-143-070
300x140x100 mm	600-143-100




Size (LxWxH)



300x140x40 mm	600-144-040
300x140x70 mm	600-144-070
300x140x100 mm	600-144-100








Size (LxW)




300x140 mm	600-145-141
	600-145-142
	600-145-143
	600-145-144
	600-145-145
	600-145-146

STERILIZATION CONTAINER

		Size (LxWxH)	
Screen Basket with Lid		271x117x45mm	600-267-050
		Lid only	271x177mm
		Bottom only	271x117x40mm
Rippled Wire Mesh		265x120x45mm	600-265-045
		Drape Retainer	280x125mm

TELESCOPE CONTAINER




STERILIZATION CONTAINER

SAFE+ CLEAN ENDOSCOPE PLASTIC CONTAINER with silicone bars on base and lid


600-505-000

Endoscopes	1 (∅ 1-10mm)
Endoscope length	max. 270mm
Dimensions	271 x 70 x 35mm




600-505-002

Endoscopes	max. 2
Endoscope length	max. 190mm
Dimensions	271 x 70 x 35mm



600-506-002*

Endoscopes	max. 2
Endoscope length	max. 390mm
Dimensions	394 x 70 x 35mm




600-501-000

Endoscopes 1

600-501-002

Endoscopes max. 2


Endoscope length	max. 450mm
Dimensions	457 x 70 x 35mm



ENDOSCOPE & INSTRUMENT CONTAINERS


600-508-000 for scopes and light guide cables

Endoscopes	max. 3 (Ø 5-10mm)
Fiber Optic Light Cable	1
Endoscope length	max. XXXmm
Dimensions	465 x 212 x 51mm



600-464-229 for fixation of instruments with rubber ties

Dimensions	464 x 229 x 51mm
------------	------------------



600-511-000 for instruments with rubber ties on 2 levels

Dimensions	502 x 239 x 102mm
------------	-------------------




STERILIZATION CONTAINER

ENDOSCOPE CLEANING BASKETS & CONTAINERS

with silicone bars on base and lid


600-292-052*

Endoscopes	max. 2
Endoscope length	max. 190mm
Dimensions	290 x 80 x 52mm




600-462-052*

Endoscopes	max. 2
Endoscope length	max. 350mm
Dimensions	460 x 80 x 52mm




600-460-055	Endoscopes 1
600-460-056	Endoscopes max. 2
Endoscope length	max. 340mm
Dimensions	455 x 86 x 43mm



600-450-070

Endoscopes	max. 2
Endoscope length	max. 330mm
Dimensions	450 x 70 x 70mm



ARTHROSCOPY INSTRUMENTS CONTAINER


600-513-000

Container for 8 instruments in upright position
with utility space for extra instruments




600-513-005

Container Case with inner tray
for 5 Chondro Picks + Hammer
or 6 Hand Grip Instruments
Dimensions: 400 x 45 x 250mm



600-513-001

Container for Arthroscopy for 16 instruments
up to 4mm in upright position with utility space
for extra instruments.
Dimensions: 464x190x70mm




STERILIZATION CONTAINER

LAPAROSCOPY INSTRUMENTS CONTAINER


600-512-000

Container for Laparoscopy for up to 12 Instruments
with maximum length of 320 mm.
Dimensions: 478x203x130mm




600-512-000

Container for Laparoscopy for up to 12 Instruments
with maximum length of 320 mm.
Dimensions: 478x203x130mm




600-513-004

Storage rack for shaft instruments
Dimensions: 220x232x155mm



600-513-003

Storage Rack for MIS instruments
Dimensions: 460x232x155mm



MAINTENANCE & ACCESSORIES




MAINTENANCE & ACCESSORIES



Identification label, 100 pieces per package

1 600-908-100 Standard

1 600-908-000 Small



Identification label

with hole for snap hook


56x17mm

unmarked	600-900-040	marked	600-900-108
	600-900-050		600-900-110
	600-900-060		600-900-120
	600-900-070		600-900-130
	600-900-080		600-900-140
	600-900-090		600-900-150
	600-900-100		600-900-160

Snap Hook

without retainer (without illustration)

600-900-170




600-900-210

Dispenser for identification labels


(8 small compartments and 1 double compartment)

as table or wall model



600-909-000

Plastic security seal, 100 pieces per package



520x230mm

600-281-004

440x230mm

600-281-003

380x230mm

600-281-002


220x230mm

600-281-001


Silicone Mat green,
cuttable to desired size




600-290-610
Paper Filter STANDARD, 100 pieces per package




600-480-100
Textile Filter STANDARD, 1 piece per package



600-380-100
PTFE (TEFLON) Filter, round, STANDARD



1 600-290-620
Disposable Paper Filters, 235x118mm with
Steam Indicator, 100 pcs.
2 600-480-200
Textile Filter, 210x95mm, 1 piece per package
3 600-380-200
PTFE (Teflon) Filter, 215x95mm,
1 piece per package



600-290-630
Filter Dispenser for disposable Filter,
round, 195x195mm

MAINTENANCE & ACCESSORIES



Filter Tray for all standard containers

Ø 185mm

600-290-640




Drape, green ●

600 x 400mm	600-900-301
850 x 850mm	600-900-302
900 x 600mm	600-900-309
900 x 800mm	600-900-310
1300 x 450mm	600-900-311
1300 x 600mm	600-900-306
1400 x 1000mm	600-900-303
2750 x 1250mm	600-900-304
2450 x 2450mm	600-900-305
1850 x 450mm	600-900-307
1850 x 600mm	600-900-308




Measuring Jug with feet, graduated

capacity	height	
0.5 l	110mm	600-098-205
1.0 l	140mm	600-098-100
1.5 l	160mm	600-098-150
2.0 l	200mm	600-098-220



Needle Sterilizing Case
Ø 65mm

Complete case	600-100-065
Perforated Interior Box only	600-100-070
Non Perforated Exterior Box only	600-100-071



SAUERBRUCH

Needle sterilising case, 12 compartments
Lid / Bottom perforation

210 x 125 x 17mm 600-210-125



Needle sterilising case, 7 compartments,
rectangular, Bottom perforation

150 x 95 x 13mm 600-115-100



Needle sterilising case, 3 compartments,
rectangular, Lid / Bottom perforation

150 x 95 x 13mm 600-115-000



Needle sterilizing case, with holes

65 x 42 x 8mm 600-110-065

MAINTENANCE & ACCESSORIES



Laboratory Bowl

∅	capacity	height	
140mm	0.5 l	55mm	600-170-151
150mm	0.75 l	75mm	600-170-165
160mm	1.0 l	75mm	600-170-175
220mm	2.0 l	75mm	600-170-220
240mm	3.0 l	95mm	600-170-240
260mm	4.5 l	120mm	600-170-275
300mm	6.0 l	120mm	600-170-310
320mm	8.0 l	140mm	600-170-325
380mm	11.0 l	130mm	600-170-380
400mm	14.0 l	155mm	600-170-415



Stainless Steel Bowl,
round

∅	capacity	height	
40mm	0.02 l	19mm	600-231-040
60mm	0.06 l	30mm	600-231-060
85mm	0.16 l	40mm	600-231-080
110mm	0.40 l	55mm	600-231-100
150mm	1.00 l	75mm	600-231-150



Instrument Box with lid

LxWxH	
100x60x40mm	600-310-100
120x80x40mm	600-310-120
160x100x30mm	600-310-160
180x120x50mm	600-310-180
200x100x50mm	600-310-200
220x150x50mm	600-312-220
260x150x50mm	600-310-260
300x200x50mm	600-310-300
400x200x80mm	600-312-400
500x270x70mm	600-312-500



Dressing Tray

LxWxH
235x190x45mm

600-340-230



Instrument Tray

LxWxH	
215x155x18mm	600-113-100
295x195x20mm	600-113-110
360x240x10mm	600-113-120
390x270x10mm	600-113-130
490x320x14mm	600-113-140
445x350x12mm	600-113-150
445x325x9mm	600-113-160



Cylinder Jar

∅	height	
30mm	95mm	600-142-001
50mm	105mm	600-142-002
50mm	135mm	600-142-003
50mm	180mm	600-142-004
80mm	135mm	600-142-006
80mm	175mm	600-142-007
95mm	205mm	600-142-009

MAINTENANCE & ACCESSORIES



Cotton Wool Dispenser

Ø	capacity	height	
60mm	0.15 l	60mm	600-397-010
85mm	0.25 l	50mm	600-397-020
100mm	0.35 l	50mm	600-397-060
85mm	0.4 l	85mm	600-397-080
85mm	0.6 l	125mm	600-397-030
100mm	0.75 l	100mm	600-397-040
120mm	1.25 l	120mm	600-398-120
150mm	1.75 l	100mm	600-398-130
150mm	2.5 l	150mm	600-398-150
200mm	3.0 l	100mm	600-398-160
200mm	6.0 l	200mm	600-398-170



Sputum Cup with lid

Ø	capacity	height	
75mm	0.3 l	75mm	600-662-075
80mm	0.4 l	85mm	600-662-080



Gynecological Bedpan
350 x 290 x 90mm
(L x W x H)

Bedpan	600-665-001
Bedpan with side handles	600-665-002
lid only	600-665-003



Irrigator Bowl
with tube fixing

∅	capacity	height	
100mm	1.0 l	165mm	600-676-010
120mm	1.5 l	175mm	600-676-015
120mm	2.0 l	195mm	600-676-020



Kidney Bowl


Length	capacity	
170mm	0.25 l	600-673-170
205mm	0.35 l	600-673-200
250mm	0.50 l	600-673-250
275mm	0.75 l	600-673-275



Slipper Pan with handle
LxWxH
360x255x90mm

Slipper Pan	600-740-300
lid only	600-740-301


MAINTENANCE & ACCESSORIES



STRINGER

Instrument Holder

102 x 70mm (4")	660-160-100
152 x 70mm (6")	660-160-150
203 x 70mm (8")	660-160-200
254 x 70mm (10")	660-160-250
305 x 70mm (12")	660-160-300
356 x 70mm (14")	660-160-350
406 x 70mm (16")	660-160-400
457 x 70mm (18")	660-160-450
508 x 70mm (20")	660-160-500



WEINSTEIN

Instrument Holder

63 x 152mm (6")	660-163-150
63 x 200mm (8")	660-163-200
63 x 250mm (10")	660-163-250
63 x 300mm (12")	660-163-300
63 x 350mm (14")	660-163-350
63 x 400mm (16")	660-163-400
63 x 450mm (18")	660-163-450



MAJO

Safety Pin, 140mm, 5¹/₂", with ball

without ball	660-100-100
with ball	660-100-140



BUNT

Safety Pin

120mm, 4 ³ / ₄ "	660-102-120
--	-------------

**ROGGE**

Sterilizing Forceps

235mm, 9¹/₄"

600-087-230

**CHEATLE**

Sterilizing Forceps

270mm, 10²/₃"

600-095-270

**DAVIS**

Sterilizing Forceps

200mm, 8"

600-080-200



Cleaning Brush for Ø 5mm

500mm, 19³/₄"

300-135-905



Cleaning Brush for Ø 10mm

500mm, 19³/₄"

300-135-910

MAINTENANCE & ACCESSORIES



Cleaning Brush

Ø3mm, 245 x 65mm	300-245-000
Ø3mm, 420 x 65mm	300-420-065
Ø7mm, 285 x 75mm	300-285-000
Ø8mm, 270 x 65mm	300-270-000




Cleaning Brush

Ø12mm, 270 x 65mm	300-215-050
Ø13mm, 215 x 60mm	300-215-000
Ø15mm, 215 x 70mm	300-215-070
Ø20mm, 215 x 65mm	300-215-080
Ø25mm, 215 x 65mm	300-215-090
Ø30mm, 215 x 60mm	300-215-100



Cleaning Brush

Ø12mm, 275 x 65mm	600-091-012
Ø15mm, 275 x 80mm	600-091-015
Ø30mm, 300 x 120mm	600-091-030
Ø50mm, 320 x 130mm	600-091-050



Cleaning Brush

Ø80mm,
490 x 180mm

300-092-080



Cleaning Brush

Ø50mm,
330 x 130mm

600-092-050



Cleaning Brush

Ø8/20mm,
300 x 190mm

300-215-110



Cleaning Brush

Ø18/33mm,
300 x 110mm

300-215-111

MAINTENANCE & ACCESSORIES



Cleaning Brush

Ø15/25mm,
320 x 70mm

300-215-112



Cleaning Brush

Ø15/25mm,
320 x 70mm

300-215-112



RZ Brush Dispenser
for 20 brushes with suspension hook


240 x 275 x 70mm

600-950-100



Cleaning Brush with plastic handle
10 pcs. / package

600-050-000



RZ Aquajet Water Pistol, complete
9 tips, 1m pressure tubing, strap,
suspension hook

600-100-000



Oil Dispenser

140mm, 5¹/₂"

700-900-250



Special Cotton for polishing metal sheaths
and instruments

700-900-220

MAINTENANCE & ACCESSORIES



Oil Spray

400ml


700-900-149



Instrument Oil

50ml

700-900-160



Tape

green 900-950-001

yellow 900-950-002

blue 900-950-003

red 900-950-004

brown 900-950-005

white 900-950-006

orange 900-950-007

black 900-950-008

purple 900-950-009



RZ Medizintechnik GmbH

Unter Hasslen 20
78532 Tuttlingen, Germany

Phone: +49 7462 9470-0

sales@rz-medizintechnik.com
www.rz-medizintechnik.com

